

Super Safety Silk Tracheal Tube

Tracheal Tube


Made of NON-DEHP PVC with high volume, low pressure cuff

Nasal/oral

Sizes: I.D. 3.0–10.0 mm

Tissue-friendly thermosensitive material, especially advantageous for long-term intubation, soft smooth tip, continuous X-ray marker, bonded connector, low pressure sealing, Safety Silk surface, black ring for intubation depth visualisation.

- Not made with natural rubber latex
- Sterile
- Single use
- Not made with DEHP



Super Safety Silk Tracheal Tube

RÜSCH

REF.	ORDER SIZE/I.D.	O.D.	CUFF Ø	LENGTH	QTY
112682 Murphy, 1 eye	3.0 mm	5.3 mm	8.0 mm	170 mm	10
	3.5 mm	6.0 mm	10.5 mm	190 mm	
	4.0 mm	6.3 mm	10.5 mm	220 mm	
	4.5 mm	6.7 mm	13.0 mm	230 mm	
112680 Magill	5.0 mm	7.3 mm	16.5 mm	250 mm	
	5.5 mm	8.0 mm	18.5 mm	280 mm	
	6.0 mm	8.7 mm	20.5 mm	290 mm	
	6.5 mm	9.3 mm	24.0 mm	300 mm	
	7.0 mm	10.0 mm	26.0 mm	320 mm	
	7.5 mm	10.7 mm	26.0 mm	330 mm	
	8.0 mm	11.6 mm	28.0 mm	340 mm	
	8.5 mm	12.4 mm	28.0 mm	345 mm	
	9.0 mm	13.0 mm	28.0 mm	350 mm	
	9.5 mm	14.2 mm	29.0 mm	350 mm	
	10.0 mm	15.0 mm	29.0 mm	350 mm	

Pedi Safety Silk Tracheal Tube

Clear Tracheal Tube

Made of NON-DEHP soft PVC

Part of the Safety Silk family

Nasal/oral (Magill)

Sizes: I.D. 2.0–5.5 mm

Especially advantageous for long-term intubation, semi-seated connector, cupped smooth black tip, continuous X-ray marker, soft version, graduated

- Not made with natural rubber latex
- Sterile
- Single use
- Not made with DEHP



Pedi Safety Silk Tracheal Tube

RÜSCH

REF.	ORDER SIZE/I.D.	O.D.	LENGTH*	BLACK TIP LENGTH	QTY
100480	2.0 mm	3.3 mm	140 mm	20 mm	10
	2.5 mm	4.0 mm	160 mm	20 mm	
	3.0 mm	4.7 mm	180 mm	20 mm	
	3.5 mm	5.3 mm	200 mm	20 mm	
	4.0 mm	6.0 mm	220 mm	30 mm	
	4.5 mm	6.7 mm	240 mm	30 mm	
	5.0 mm	7.3 mm	260 mm	40 mm	
	5.5 mm	8.0 mm	280 mm	40 mm	

*Length without connector