

## Spiral-Flex Reinforced Tracheal Tube

High-volume, tapered,  
low-pressure cuff

Nasal/oral, with cuffed version

Sizes: I.D. 5.0–9.0 mm

Reduced risk of kinking due to stainless steel, spiral wound reinforcing wire within the tube wall; 15 mm connector is permanently bonded within the tube, blue cuff

- Not made with natural rubber latex
- Sterile
- Single use

### Cuffed Oral

Cuffed, oral precut length with reinforcing sleeve to minimise the risk of patient biting and occluding the tube; the oral length tube has a reinserted Sher-i-slip stylet.



### Spiral-Flex Reinforced Tracheal Tube

**SHERIDAN**

REF.	ORDER SIZE/I.D.	O.D.	CUFF Ø	LENGTH	QTY
Cuffed oral					5
5-12510	5.0 mm	7.5 mm	19.0 mm	216 mm	
5-12511	5.5 mm	8.2 mm	22.0 mm	221 mm	
5-12512	6.0 mm	8.8 mm	24.0 mm	240 mm	
5-12513	6.5 mm	9.4 mm	26.0 mm	260 mm	
5-12514	7.0 mm	10.0 mm	27.0 mm	279 mm	
5-12515	7.5 mm	10.6 mm	28.0 mm	295 mm	
5-12516	8.0 mm	11.4 mm	28.0 mm	300 mm	
5-12517	8.5 mm	12.0 mm	29.0 mm	321 mm	
5-12518	9.0 mm	12.6 mm	29.0 mm	321 mm	
Cuffed oral/nasal					5
5-12610	5.0 mm	7.5 mm	19.0 mm	246 mm	
5-12611	5.5 mm	8.2 mm	22.0 mm	276 mm	
5-12612	6.0 mm	8.8 mm	24.0 mm	286 mm	
5-12613	6.5 mm	9.4 mm	26.0 mm	297 mm	
5-12614	7.0 mm	10.0 mm	27.0 mm	306 mm	
5-12615	7.5 mm	10.6 mm	28.0 mm	318 mm	
5-12616	8.0 mm	11.4 mm	28.0 mm	327 mm	
5-12617	8.5 mm	12.0 mm	29.0 mm	327 mm	
5-12618	9.0 mm	12.6 mm	29.0 mm	327 mm	