


LMA Protector Airway

Laryngeal Mask with dual gastric access

Made of silicone with Cuff Pilot Technology option

The LMA Protector Airway is a single use second generation device which features an innovative dual gastric access channel, and a true Second Seal Technology with the upper oesophageal sphincter. It allows direct intubation using visual guidance.

- Not made with natural rubber latex
- Sterile
- Single use



LMA Protector Airway

LMA

REF. WITH CUFF PILOT TECHNOLOGY	REF. WITH PILOT BALLOON	MASK SIZE	PATIENT WEIGHT	MAXIMUM INTRACUFF PRESSURE**	MAXIMUM ETT ID (MM)	LARGEST SIZE OG TUBE	QTY
192030	195030	3	30–50 kg	60 cm H ₂ O	6.5	16 Fr.	10
192040	195040	4	50–70 kg	60 cm H ₂ O	7.5	18 Fr.	10
192050	195050	5	70–100 kg	60 cm H ₂ O	7.5	18 Fr.	10

ETT=ENDOTRACHEAL TUBE | OG=OROGASTRIC TUBE

* LMA Protector Airway with Cuff Pilot Technology only.

** It is recommended that the cuff be inflated to a maximum intracuff pressure of 60 cm H₂O.