

Laser Tube

Preformed Tracheal Tube

With defined laser-resistance made of white soft rubber

For laser intervention in the laryngeal space

Oral (Magill), 1 eye (Murphy)
 Sizes: I.D. 4.0–8.0 mm

Fixed white connector, two-way stop-cocks with Luer and Luer-lock adapter and locking mechanism, cupped smooth tip, two cuffs inside one another, two pilot balloons, non-kinking shaft, Laserguard film 17 cm long consisting of Merocel foam and corrugated silver foil

- Sterile
- Single use



Laser Tube						RÜSCH
REF.	ORDER SIZE/I.D.	O.D.	CUFF Ø OUTSIDE	CUFF Ø INSIDE	LENGTH*	QTY
102004	4.0 mm	8.0 mm	20 mm	20 mm	400 mm	2
	5.0 mm	9.0 mm	20 mm	20 mm	400 mm	
	6.0 mm	10.8 mm	24 mm	24 mm	400 mm	
	7.0 mm	12.3 mm	28 mm	28 mm	400 mm	
	8.0 mm	13.4 mm	28 mm	28 mm	400 mm	

*Length without connector



Caution:
 These products contain natural rubber latex which may cause allergic reactions.