



*EZ-Blocker Endobronchial Blocker
combined with ETT
(not included in set)*

EZ-Blocker Set

Y-Shaped Endobronchial Blocker

For lung isolation and one-lung ventilation

- For surgery with fiberoptic or video bronchoscopes
- 75 cm long, exact positioning via an endotracheal tube
- Bifurcated distal extensions, both ends with an inflatable cuff and a central lumen, colour coded, minimal risk of dislocation during procedure
- Radiopaque shaft with depth markers, cuffs made of polyurethane.
- EZ-Multiport Adaptor to facilitate ventilation, EZ-Blocker Endobronchial Blocker placement and the introduction of fiberoptic or video bronchoscopes and suction catheters
- Not made with natural rubber latex
- Sterile
- Single use
- Not made with DEHP



EZ-Blocker Set

RÜSCH

| REF. | SET COMPONENTS | LENGTH | SIZE | QTY |
|--------------|--|--------|------|-----|
| MG-02770-002 | 1 EZ-Blocker Endobronchial Blocker 1 EZ-Multiport Adaptor with scope lid 1 Oxygen adaptor 2 Dust caps | 75 cm | 7 Fr | 1 |

