

A-7000 / A-8000 SERIES


SAFE AND ACCURATE

A-7000 / A-8000 – WET SUCTION /WET SEAL

The PLEUR-EVAC wet suction series has been the industry standard since its introduction in 1967. Wet suction models are safe and accurate. The A-7000 series offers a patented patient air leak meter to quantify the patient air leak, needleless tubing connectors for collection sampling and is available with premated or optional ATS bags.

- **WATER-CONTROLLED SUCTION COLUMNS**
provide safe suction levels
- **HIGH NEGATIVITY FLOAT VALVE**
vents loss of water seal
- **HIGH NEGATIVITY RELIEF VALVE**
permits manual filtered venting of excessive negative pressure

WET SUCTION / WET SEAL – ALL MODELS LISTED BELOW CONTAIN ATS CONNECTORS AND A PATENTED AIR LEAK METER



A-7000 SERIES – ADULTS/CHILDREN

PLEUR-EVAC

REF.	DESCRIPTION	VOL. COLLECTION CHAMBER	QTY
A-7000-08LF	single collection chamber	2500 ml	6
A-7002-08LF	dual collection chamber	950 and 1900 ml	6
A-7050-08LF	with premated ATS Bag	2500 ml	3

A-7000 ATS BAG

PLEUR-EVAC

REF.	DESCRIPTION	QTY
A-1500-08LF	PLEUR-EVAC ATS Bag, fits all 7000 series models	6



S-1500-08LF

WET SUCTION / WET SEAL – THESE SERIES DO NOT HAVE ATS CAPABILITY

A-8000 SERIES – ADULTS/CHILDREN

PLEUR-EVAC

REF.	DESCRIPTION	VOL. COLLECTION CHAMBER	QTY
A-8000-08LF	single collection chamber	2500 ml	6
A-8002-08LF	dual collection chamber	950 and 1900 ml	6



A-8000-08LF

A-8000 SERIES – INFANTS

PLEUR-EVAC

REF.	DESCRIPTION	VOL. COLLECTION CHAMBER	QTY
A-8020-08LF	Infant Model	150 ml	3