

Accessories for Serres suction system

- All accessories are compatible with Serres suction bags and are easy and fast to connect.
- **Collection cup** is designed for collecting particles from the suctioned liquid.
- **Measuring cup** for accurate measuring of small liquid volumes.
- **Splash-vac** is ideal for removing large amount of liquids from the floor conveniently.
- **Bacterial filter** can be connected between suction canister and suction pump.
- **Smoke filter** prevents premature filter blockage by smoke.
- **Vacuum source protector** prevents the contamination of the hospital's central vacuum system and suction pump in case of an unpredictable failure of the fluid collection system.

Product no.	Product name	Packing size (pcs)
57247	Collection cup 300ml	108 (18x6)
57242	Collection cup 300ml, without angle connector	108 (18x6)
57257	Measuring cup 250ml	72 (12x6)
57270	Splash-vac + 3m suction tubing	25
57702	Bacterial/Viral filter	10
57710	Smoke filter	10
6001114	Vacuum source protector	15



57247



57257



57270



57702



57710



6001114

Serres transport box

Serres suction bags are foldable to fit up to four times more bags into a box compared to similar products on the market, leading to more efficient transport and space savings in storage.



Each transport box holds 6 dispenser boxes

Each dispenser box holds 24 standard 3L Serres suction bags

One transport box fits 216 one-liter or 144 two- or three-liter Serres suction bags in total