

FLUID COLLECTION

Valves and connectors for Serres canisters

- Serres valve is used to open and close suction when the smoke filter is used.
- On-off valve can be used to open and close suction. It replaces the grey angle connector.
- T-connector is used in serial connections and can replace the Serres valve when the smoke filter is used
- The reusable Angle connector is available as a spare part if it is disposed of accidentally.

Product no.	Product name	Packing size (pcs)
57958	Serres valve	10
57959	On-Off valve	10
57610	T-connector	100
57618	Angle connector	20



57958

57959



57610

57618

Solidifying agent

The Serres solidifying agent ensures safe handling of fluid waste, minimizing risk of spillage or contact with potentially infectious fluid waste.

Pre-gelled Serres suction bags

Pre-gelled Serres suction bags with solidifier enable safer handling.

Product no.	Product name	Packing size (pcs)
57557	Pre-gelled bag, 1000ml	192 (32 x 6)
57567-0	Pre-gelled bag, folded, 2000ml	132 (22 x 6)
57587-0	Pre-gelled bag, folded, 3000ml	120 (20 x 6)



Solidifying agent

Extra solidifier is available in separate pouches or in a 2kg tub. Solidifier can be added into the suction bag through a serial port before or after the procedure.

Product no.	Product name	Packing size (pcs)
57550	Solidifying agent, 25g	24/288 *
57551	Solidifying agent, 15g	20/200 *
57553	Solidifying agent tub, 2kg	1

* Pouches are packed in plastic bag of 20/24 pcs. The total number of pouches in a storage box is 200/288 pcs

