

Suction canisters

Reusable Serres canisters need replacement very rarely, as they are extremely durable. Serres suction canisters are certified as measuring devices with an accuracy of +/- 100ml. Canisters are equipped with built-in brackets for mounting on walls, rail supports or trolleys. Canisters include reusable angle connectors for vacuum tubing. The compact oval-shaped 1000ml canisters are ideal for tight places in ambulances, intensive care units and anaesthesia work stations.

Suction canisters, transparent

Transparent canister with wide scale allow for a clear view of suctioned liquid.

Product no.	Product name	Packing size (pcs)
57308	Suction canister, 1000ml	45
57328	Suction canister, 2000ml	24
57348	Suction canister, 3000ml	16



Suction canisters, blue

Blue canisters are designed for sensitive circumstances to ensure comfort of patients and families.

Product no.	Product name	Packing size (pcs)
57309	Suction canister blue, 1000ml	45
57329	Suction canister blue, 2000ml	24



Flushing canister

In intermittent use, suction tubes can be cleaned and kept operational by suctioning small amount of water from the Flushing Canister.

Product no.	Product name	Packing size (pcs)
57319	Flushing canister, white, 1000ml	45

